
*Flat  
Adenoma  
Resection  
Instrument  
(TYP C)*

# Flat Adenoma Resection Instrument [Typ C]

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The flat adenoma resection instrument (type C) was developed to resect especially large polyp-like and flat lesions > 2 cm to approx. 4 cm en bloc close to the muscularis propria and with low mechanical and / or thermal artefacts in a relatively short time (pure resection time, i.e. without the time for marking, undercutting and resection). (pure resection time, i.e. without time for marking, undercutting and resection < 1 minute).

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
## The essential features

The objective of developing new procedures and instruments for endoscopic removal of large polyps (in particular, those exceeding >2 cm in diameter) and flat lesions of the mucosa of the gastrointestinal tract is to fulfil the relevant requirements of oncology, namely the removal of pathological tissue. This is to be achieved by removing pathological tissue in sano and pathology, namely the removal of pathological tissue, preferably en bloc and below SM1, or close to the muscularis propria, optimally according to the incision A (see picture).

The currently available methods of endoscopic polypectomy (ESR) and mucosal resection (EMR), which are performed with conventional polypectomy and HF surgical resection loops, respectively, are unable to reliably fulfil the aforementioned requirements of oncology and pathology.

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A procedure called Endoscopic Submucosa Resection (ESR) and appropriate instruments (Flat Adenoma Resection Instrument) are available for clinical application after extensive in vitro testing.


The most important feature of these instruments is the effector, which is available as a symmetrically opening and closing HF surgical resection snare.

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
### Symmetric RF-resection snare.

- 1 · catheter
- 2,3 · electrically insulated parts
- 4 · RF-cutting wire
- 5 · electrically insulated spur
- 6 · emitransparent distal end of the catheter
- 7,8 · markings which can be seen endoscopically through the distal end 6 of the catheter for controlling the movement of the snare especially during the RF-cutting.

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### Results of an ESR with a Flat Adenoma Resection Instrument

(Resector: OA Dr. med. Dr. Gölder (III. med. clinic of the Central Hospital Augsburg, CA Prof. Dr. H. Messmann) confirms that the above-mentioned requirements of oncology and pathology have been fulfilled.

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## Technical Data

Opening and closing of the effector	symmetrical
Rotating the snare	at the catheter axis
Form des Effektors	see picture
Maximum opening width of the snare	3 cm
effektive Länge des Schneidrahts	15 mm
Diameter of the cutting wire	0,39 mm
Length of the catheter	200 cm
External diameter of the catheter	2,3 mm

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Feinwerktechnik GmbH