



Loops



S-L(D)Ø-a x b -sym



S-L(D)Ø-a x b -hex

Construction

Loop
bent, lasered or brazed strand made of stainless steel.

Pull pope
strand of stainless steel or rather wire with pull rod made of stainless steel.

Explanations
 LØ- diameter of strand
 DØ- diameter of wire
 a- width
 b- length
 sym- symmetric
 hex- hexagonal

Customer requirements
 Besides our standards, we would be pleased to be your manufacturer when it comes to individual specifications and wishes for Loops.